



PRODUCT DATASHEET



High Speed Digital EMCCD Camera for Bio-Imaging

Introducing the Rolera EM-C²
EMCCD camera – a versatile
bio-imaging platform which combines
unprecedented 40 MHz speed with
high resolution and low light sensitivity.

The EM-C² boasts a brand new feature called Easy-EM™, the one-click optimization of your camera's EM gain settings in all modes of speed. Instant experimental repeatability is possible with Easy-EM, as each camera is uniquely characterized and operational with the best read noise reduction, dynamic range, and limited EM noise injection.

Stream video of your experiment at 34.2 full frames per second. View live cells in high speed, with 200+ FPS using binning and ROIs.

Get up and running fast with intuitive hardware triggering options, compatibility in popular microscopy control and imaging software suites.

Speed is guaranteed with the Rolera EM-C²'s dedicated IEEE 1394b

FireWire interface.

High Speed, High Sensitivity Digital EMCCD Camera



| teatures | penetits | | |
|------------------------------------|--|--|--|
| Introducing Easy-EM [™] | Each camera is individually characterized to determine its optimal EM Gain setting. Easy-EM sets your camera to this highest-sensitivity mode. | | |
| Active cooling | Actively cooled to -50° to reduce thermal noise | | |
| High-Speed Readout | Previewing and focusing in real time 165.8fps with 8x8 binning 34.2fps full resolution @ 14 bits Ideal for automated imaging applications | | |
| Low-Noise Electronics | Quantitation and imaging of low light levels | | |
| Flexible Exposure Control | Optimal integration over a wide range of light levels | | |
| External Sync and Trigger | Tight synchronization with flashlamps, automated filters, shutters, and microscope stages | | |
| Binning | Increases sensitivity for quantitation and imaging of very low light levels Increases frame rate | | |
| IEEE-1394b FireWire 800 Connection | Simple connectivityEase of use and installation | | |
| 14-bit readout | ■ More than 16,000 gray levels per pixel | | |
| Integrated digitization | ■ No RF interference | | |
| No need for external controller | ■ Keeps bench space clutter-free | | |
| Easy installation | ■ Ready to use in less than 10 minutes | | |
| Software included | ■ Enables a variety of life-science applications | | |
| RGB module available | ■ High-resolution color images | | |
| applications | | | |

applications

The Rolera EM-C² is the fastest, highest resolution EMCCD camera which performs in a variety of low-light bio-imaging applications.

- Live cell imaging
- Spinning disk confocal microscopy
- Fluorescent protein imaging (BFP, GFP, YFP, RFP, FRAP)
- TIRF
- Ratiometric Imaging
- Ion Imaging

Rolera EM-C² Specifications

| ccd sensor | | | |
|---|--|--|--|
| Light-Sensitive Pixels | 1004 x 1002 | | |
| Binning Modes | 2, 4, 8 | | |
| Exposure/Integration Control | 200μs to 17.9min | | |
| Sensor Type | Texas Instruments TX285 Frame Transfer EMCCD | | |
| Pixel Size | 8µm x 8µm | | |
| Linear Full Well | 32,700e- | | |
| Read Noise (Typical) | 20e- rms @ 40MHz 20e- rms @ 20MHz 18e- rms @ 10MHz Read noise effectively reduced to <1 e- rms with EM gain enabled | | |
| Dark Current | 0.06 e-/pix/s | | |
| Cooling | -50°C (regulated) | | |
| Digital Output | 8 bits/14 bits | | |
| Readout Frequency | 40, 20, 10MHz | | |
| Frame Rate | 34.2fps full resolution @ 14 bits (40MHz) | | |
| camera | | | |
| Computer Platforms/ Operating Systems* | Windows® 7, Windows XP, Windows Vista (64/32 bit) | | |
| Digital Interface | IEEE-1394b FireWire | | |
| External Trigger | TTL Input | | |
| Trigger Types | Internal, Software, External (Edge-Hi/Edge-Low/Pulse-Hi/Pulse-Low/Strobe-Hi/Strobe-Low) | | |
| External Sync | TTL Output | | |
| External RGB Filter Control | Support for RGB filter | | |
| Optical Interface | 2/3", C-mount optical format | | |
| Threadmount | 1/4" – 20 mount | | |
| Power Requirements | 40 watts at 12 volts | | |
| Weight | 952g (2.1lbs) | | |
| Warranty | 2 years | | |
| Operating Environment | 0 to 25°C, 80% relative humidity non-condensing | | |
| Storage Temperature | | | |

cell imaging frame rates

No BINNING providing full spatial resolution

| | Approx. Size* (um) | Required Region of Interest | Frames/Sec | |
|-------------------|--------------------|-----------------------------|------------|--|
| Endosome/Lysosome | 5x5 | 38x38 | 205 | |
| Cell Nucleus | 10x10 | 75x75 | 171 | |
| Golgi | 15x15 | 113x113 | 147 | |
| Small Cell Region | 20x20 | 150x150 | 134 | |
| СНО | 25x25 | 188x188 | 124 | |
| HeLa | 30x30 | 225x225 | 111 | |

^{*}Using average cell sizes with a 60X objective to calculate ROI requirements

included

- Rolera EM-C² Bio-Imaging Microscopy Camera Model: 01-ROL-EMC2-R-F-M-14-C (monochrome)
- Power supply
- *IEEE1394b FireWire cable* (9-pin to 9-pin)
- IEEE1394b FireWire PCI Express Card
- QCapture Suite software for PC
- Limited Warranty

camera options

- Performance Assurance Program (extended warranty)
- RGB Filter Module

spectral response









